



1. Proper care and periodic maintenance to your instruments will provide long-term use and enable you to perform the best possible patient care.
2. The user is responsible for inspecting the instruments prior to use.
3. Daily instruments inventory check will help to track any damages or lose.
4. The instruments should be handled with care.
5. Store instruments in a dry and dust free place after sterilization.
6. Using instruments for purposes other than what they are designed for leads to damages and breakage.
7. Defected Instruments should be taken out of service immediately, for example; instruments that are dull, shows any sign of corrosions, misaligned, defected in any form or instruments with weakened blades.
8. Always Keep your hygiene scalars and curettes sharp all the time in order to facilitate the removal of debris and avoid excessive pressure on the instrument tip that may lead to breakage.
9. Engravings added to the instruments surface by the user is discouraged as it is considered a potential inroad for oxidization and/or corrosion.
10. Return any damaged instruments to the Dental Instruments Coordination Office for exchange.

CONTACT US :

Tel: +965-25627862

Fax: +965-25613845

e-mail: dent.instruments@hotmail.com

twitter: [@dental_store](https://twitter.com/dental_store)